realcheck™ PRO 2

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1. INTRODUCTION

1.1 Features of the REALCHECK PRO 2

REALCHECK PRO 2 features QMS technology (Quick Measuring System) for a quicker and easier reading of pressure parameters.

QMS - "Initially the cuff inflates to detect the right pressure on the arm, then cuff inflation stabilizes and goes on steadily. After detecting the blood pressure value the unit reads the Diastolic and Systolic values. As soon as the Systolic value has been measured, the inflation pump stops deflating the cuff rapidly, thus shortening the measuring time and improving patient comfort during measurement."

Once the reading has been taken, the rapid deflating process commences and the display shows the systolic and diastolic pressure together with the heart rates. The high precision has been clinically tested and the system has been specially designed for the utmost ease of use.

REALCHECK PRO 2 features a large LCD display that permits clear reading of the functional status and cuff pressure throughout the measurement process.

Before using, read the contents of this manual carefully and then store it in a safe place. Contact a physician for further queries concerning blood pressure and how to measure it.

Attention!

1.2 Important information about self-measurement

- Do not forget: **self-measurement means Control**, not diagnosis or treatment. Unusual values must always be discussed with the doctor. Any change in dosage of medications prescribed by the doctor should be discussed with him/her personally.
- The pulse display is **not** suitable for checking the frequency of heart pacemakers!
- In cases of cardiac irregularity (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with the doctor.

2. INFORMATION REGARDING BLOOD PRESSURE AND HOW TO MEASURE IT CORRECTLY

2.1 Hypertension

Hypertension is one of the most important vascular risk factors and is directly linked to the main causes of death in the western world (heart attack and ictus). Checking the arterial pressure must then be considered an indication of the health of the cardiovascular system. Like all biological variables, blood pressure is changeable and is influenced by numerous parameters which can be internal, (hormones, stress, anxiety etc.) or external (temperature, smoke, coffee etc.), to our organism; it indicates the pressure that the blood, pumped by the heart, exerts on the arteries: *maximum or systolic pressure – minimum or diastolic pressure.*


Fig. 1

Source: Guideline (1999) World Health Organization (WHO)
- International society of Hypertension (ISH)

It tends to increase with age and is influenced by the seasons, the climate and the altitude. By measuring it (in millimetres of Mercury – mmHg -) normal conditions can be established, as well a situation of hypertension.

Variations in blood pressure in response to diverse stimuli during the day make it hard to judge on the basis of a single measuring. It is consequently advisable to measure it twice a day (once in the morning on waking and once in the late afternoon).

Persistent hypertension should be treated with medications under the strict supervision of a doctor, and certain habits should be corrected (salt consumption, a sedentary way of life, obesity, smoke, alcohol, stress).

The measuring of blood pressure enables you to verify the effectiveness of the measures adopted to keep the situation under control. Monitoring your own blood pressure is a practical and easy-to-repeat process when reliable devices are used and has the added advantage of reducing the feeling of unease often experienced by patients having their blood pressure taken in doctors' surgeries or hospitals, which can result in altered values. However it is important to underline that the option to measure your own blood pressure should not become a reason to avoid a visit to the doctor, and it certainly should not lead to 'do-it-yourself' treatment.

It is therefore essential that frequency of measuring, the results of the measuring and eventual therapies should be agreed with your doctor. It is always unwise to alter doses, or drugs prescribed by your doctor, on your own initiative.

**READ INSTRUCTION MANUAL CAREFULLY BEFORE USE
THE DEVICE MUST NOT BE DISASSEMBLED
FOR ANY PROBLEM PLEASE CONTACT THE TECHNICAL SERVICE OF DISTRIBUTOR
FOR SPECIFIC INFORMATION ON YOUR OWN BLOOD PRESSURE PLEASE CONTACT
A DOCTOR
USE ONLY THE ORIGINAL DC ADAPTER**

3. IMPORTANT SAFETY RULES

1. Before using the device, always check the electrical requirements on manufacturer label and verify the type of battery to use itself.
2. Particular precautions must be made concerning electromagnetic compatibility. The medical device must be set-up and used according to information supplied with the present manual;
3. Respect the safety rules and particularly:
 - Children and/or disable persons may only use the device under the supervision of an adult in possession of his/her full faculties.
 - If the device is stored at low temperature, do not use immediately. Leave the device at room temperature for at least one hour.
 - Avoid extreme environmental conditions and the exposure to direct sunlight.
 - **NEVER DIP ANY PART OF THE DEVICE AND ACCESSORIES IN WATER, IT ISN'T WATER PROOF.**
 - The device must be used only by adult people.
 - Cleaning the device with soft and dry cloth. Do not use petrol, thinners or similar solvents. Spot on the cuff can be removed carefully with a damp cloth. Do not wash any part of the device.
 - If you do not use the device for a long time, please take the batteries off. **KEEP THE BATTERIES OUT OF THE REACH OF CHILDREN.**
 - It is recommended the use of alkaline batteries.
 - Do not wash or iron the cuff

- Do not squeeze the air hoses nor the cuff
 - Keep the cuff away from cutting edges and do not twist it
 - Avoid to drop the device or treat it roughly in any way.
4. The device must not be disassembled. For any problem please contact the Technical Service of the Distributor. **Failure to comply with the above invalidates the safety of the device and the guarantee claims.**
 5. **This medical device must be used exclusively for the use for which it has been designed and as described in this manual.** Any different use must be considered incorrect and therefore dangerous; the manufacturer and the distributor cannot be considered liable for damage caused to people, things or animals by improper, incorrect and/or unreasonable use.
 6. Measurement on people with serious peripheral circulation or renal problems, people who have had an ictus or if are unconscious and the case of serious arrhythmia may be difficult and bring extra error; it is recommended to consult the physician.
 7. In case of malfunctions contact the technical service of the distributor. If by technical check the device will result opened by user, the guarantee will be expired.
 8. The manufacturer and distributor cannot be considered liable for damages to people, things or animals, caused by connection to electrical installations which do not conform to safety regulations, use of adapter or extension cables not mentioned in this manual, connection to main socks with voltage and frequency not in compliance with data printed on rating label of your device, repairs or changes by user or people not authorized.

4. PUTTING THE BLOOD-PRESSURE MONITOR INTO OPERATION

4.1 Inserting the batteries

After you have unpacked your device, first insert the batteries. The battery compartment is located

on the underside of the device (see illustration).

a) Remove cover as illustrated in Fig. 2

b) Insert the batteries (4 x size AA 1,5V), **thereby observing the indicated polarity** (Fig.3) and close the cover.

Caution, the life span of the batteries included can be shorter than expected in consequence of storage period.


Fig. 2


Fig. 3

4.2 Battery replacement

1. In order to keep the record of reading and timer information, please replace batteries under the Clock mode (only date and clock on the Display not lighted).

2. If the replacements process exceeds 10 seconds, the data will be cleared out and the timer information will be reset.

3. When you install or replace the batteries over 10 seconds, the device will automatically

reset and display all symbols on the LCD for 1 seconds (Fig. 4A), then shows the information of clock (Fig. 4B). To avoid the clear of memory plug the adaptor DC during the battery replacement, and unplug it once done it


Fig. 4A


Fig. 4B

Fig. 4C


Note: If the device cannot be actuated, please remove the batteries and wait some minutes, then replace the batteries again.

4. When weak battery symbol appears on the display (Fig. 4C), please replace with four new batteries at the same time. Do not leave exhaust batteries into device compartment.
5. Do not use rechargeable batteries. Use only batteries in compliance with specification declared in this manual.
6. Remove the batteries from the compartment if the device will be used with DC adaptor for long periods.

4.3 Use with DC adaptor

This device can operate with DC adapter accessory and enclosed in the box. To operate with adaptor, you must only connect the DC connector to the plug on the right of the device and connect the adapter to a main power socket adapt to your adapter. Please use only the enclosed adapter or its original spare part. In fact only this adapter model is certified for use with you REALCHECK PRO 2. Remove the batteries from the compartment if the device will be used with DC adapter for long periods

4.4 Display Explanation

The Pulse mark ♥ also has function to advise, at the end of measurement, in case of Irregular Heartbeat (arrhythmia). When Heartbeat results irregular the Pulse mark (Fig. 5) will flash automatically at the end of measurement.

| Figure | Range | Systolic (mmHg) | Diastolic (mmHg) |
|---|--------------|-----------------|------------------|
|  | Hypertension | > 140 | > 90 |
|  | Normal-High | 130~139 | 85~89 |
|  | Normal | < 129 | < 84 |
|  | Hypotension | < 90 | < 50 |

Data in accordance with:
 (1999) WHO (World Health Organization) / ISH (International Society of Hypertension) Definitions and Classification of Blood Pressure Level
 2003 ESH(European Society of Hypertension) / ESC(European Society of Cardiology) Guidelines


*INTERPRETATION ICONS

These 4 symbols are a quick, convenient way of gauging the condition of one's blood pressure, following the table below.

4.5 Setting

Press one time the **"Mode"** key and device will enter in the setting mode for Clock, Alarm, Pc linkage settings. To any other press of **"Mode"** key you can see the different setting modes in the following sequence: Clock ad Date (with **CL** indication blinking on the display), Alarms (with **AL** indication blinking on the display), PC connection (with **PC** indication blinking on the display).


4.5.1 Clock and Date

Press the **"Mode"** key a first time to enter in the setting mode and a second time to set the Clock mode, the indication **"CL"** will blink on the display (as in the picture).

Press now the **"Set"** key to enter in the Clock mode and move along the data that you need to set and use the **"Adjust"** key to modify these.

The sequence of number appeared on the display is Year, Month, Day, Hour, Minutes. Press the **"Set"** key to the number that you want to adjust, when this number blinks on the display, you can press **"Adjust"** key to increase and change so it.



Example: after to be entered in the Clock mode as explained, if you want to adjust the **Hour**, press the **"Set"** key until when the numbers of **Hour*** blinks on the display, after press the **"Adjust"** key to increase and change it. After to have finished the Hour setting, press again the **"Set"** key to **pass to other data**.

To close the setting mode, when you have finished your setting procedures, press the **Off** key.

4.5.2 Clock Alarms

There are three sets of day alarm that you can set to remind to take medication or to measure your pressure.

key to enter in the setting mode and press it other two times to open the Alarm Mode, that will be indicated by you by **"AL"** indication blinking on the display. Now press the **"Set"** key and the device will enter in the Clock Alarm Mode and the indication


"01" will blink on the display, this means that you are in the first set of Alarm.

For the second and third set you can press the **"Adjust"** key (sequence 01 - 02 - 03 will appear blinking on the display). When one of this alarm number blink on the display you must press the **"Set"** key to select

it (when selected the number of alarm will stop to blink) and after you must press the **"Adjust"** key to active the the alarm (the symbol ☹ will appear on the display)

Now the alarms are activated and you can adjust the time (hours and minutes) using the **"Set"** key to select the data and **"Adjust"** key to modify the data.

To close the alarm settings press the **"Off"** key after to have finished your adjustments. If you

have properly activated the alarms, the symbol ☹ will be displayed. To turn off the alarm press the **"Start/Off"** key but remember that the alarms will be active for all days until you will change their setting by On to Off status, using the same procedure already described above. In any moment it is possible to modify the data memorized in the alarms, following the normal set procedure.

4.5.3 PC connection

Attention: This function is only to save and manage on your PC the measuring data recorded in the device memories. This function cannot change or manage any parameters or setting of your device.




To use correctly this function it is necessary, for the user, an experience in using of Personal Computer.

First of all, please proceed to install all software required for use and contained in the enclosed CD-rom.

Print and read very well the instructions of use and follow them carefully.

Adapt for Win98SE/ME/Win2K/XP/2003/Vista32.

To connect REALCHECK PRO 2 to your PC where you have installed the software, use the USB cable enclosed in the pack. Press the **"Mode"** key to enter in the setting mode and press again until the **"PC"** indication will blink on the display.

Now press **"Set"** key and the device will enter in the PC linkage mode and the **M1** indication will blink on the display. With **"Adjust"** key you can select the memory channel (M1 or M2) that you want download on your PC.

If the PC is connected and ready to work (as indicated in the instruction you have downloaded by CD-rom),

you can proceed with the automatic download, pressing the “Set” key. Remember that you will download the data from the Memory which blinks on the display (M1 if the “M1” mark is blinking and M2 if the “M2” mark is blinking).

5. HOW TO MEASURE YOUR BLOOD PRESSURE



5.1 Before the measurement:

- Avoid eating, smoking as well as all forms of exercise before the measurement. All these factors influence the measurement result. Try and find time to relax by sitting in an armchair in a quite atmosphere for about ten minutes before the measurement.
- Remove any garment that fits closely to your upper arm.
- Measure always on the same arm (normally left).
- Attempt to carry out the measurements regularly at the same time of day, since the blood-pressure changes during the course of the day.

5.2 Common sources of error:

Note: Comparable blood-pressure measurements always require the same conditions! These are normally always quiet conditions.

- All efforts by the patient to support the arm can increase the blood-pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower (or higher) than the heart, an erroneously higher (or lower) blood-pressure will be measured!


- Cuff which is too narrow or too short result in false measurement values. Selecting the correct cuff is very important. The cuff size depends upon the circumference of the arm (measured in the centre). The cuff range is printed on the cuff. If this is not suitable for your use, please contact the Technical Service of the distributor.

Warning: Do not use cuffs other than the original cuff contained in this pack.


- A loose cuff or a kinked tube may cause false measurement values.
- Please avoid strong electromagnetism interference, (such as mobile telephone, microwave oven, etc). This can cause low accuracy in the measurement.
- Use only alkaline batteries.
- Do not judge the measuring results and adjust your treatment yourself, it may be dangerous. Please give measuring results to physician to make a diagnose. Please take your medications as by physician's instructions.
- This equipment may stop working if there is a strong interference of electrostatic discharge.
- If the device is stored at low temperature, do not use it immediately. Leave the device at room temperature for at least one hour before using it.
- The Arrhythmia may cause mistakes in the measurements, so is good to use with great caution the device on peoples that have arrhythmia symptoms.

5.3 Fitting the cuff

a) Close the cuff with the Velcro closer in such a way that it lies comfortably and is not too tight. Connect the armband tube to the appropriate opening on the left of the apparatus as shown in the illustration.


b) Pass the end of the cuff (with sewn-in rubber stopper) through the metal stirrup so that a loop is formed. The Velcro closer must be facing outwards. (Ignore this step if the cuff has already been prepared.)

c) Push the cuff over the left upper arm so that the tube points in the direction of the lower arm.


d) Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately 2 to 3 cm above the elbow and the air tube lays on the inner side of the arm.

Important! Ensure that the reference sign (symbol Φ) is over the Brachial Artery that is placed on the inside of the arm.

e) Tighten the free end of the cuff and close the cuff with the closer.

- f) There must be no free space between the arm and the cuff as this would influence the result. Clothing must not restrict the arm. Any piece of clothing which does (e.g. a pullover) must be taken off.
- g) Lay the arm on the table (palm upwards) so that the cuff is at the same height as the heart. Make sure that the tube is not kinked.


- h) Remain seated quietly for two minutes before starting the measurement.
- i) If the "INDEX" line falls within the OK range, indicated at the edge of the cuff, this cuff is suitable for the use. If the "INDEX" line falls out of the OK range, it could be necessary a special cuff.

Note:


If it is not possible to fit the cuff to the left arm, it can also be placed on the right one. However all measurements should be made using the same arm.

5.4 Measuring procedure

1. Press **Start** button to operate the device.
2. The display will switch on with the indication 00 blinking.
3. About 2-3 seconds later, the pressure will be reset, then the cuff begins to inflate and the pressure value will be displayed.

4. This monitor will inflate automatically if the system detects that your body requires more pressure for measurement.
5. The cuff pressure will fall slowly during measuring. Do not shake the cuff or the monitor. Keep quiet.
6. Along with the measurement, besides the values of blood pressure which can detect if your blood pressure has problem or not, the heartbeat rate is also recorded (along with the measurements, the heart sign will appear on the display, blinking to any heartbeat) and detected if you have an irregular heartbeat(arrhythmia) or not, the heartbeat reading is indicated on the display at the end of measurement.
7. After finishing measurement the cuff deflates rapidly and automatically and the values of systolic and diastolic pressures and pulse rate will be displayed on the display together an automatic interpretation symbol (as indicated at paragraph 4.4 of this instruction manual) and the heart symbol (blinking or fix as below pictures).

Note: It is recommended to consult a physician if one of following symbols appear blinking on the display at the end of measurement:


Hypertensive values was detected

Hypotensive values was detective

An irregular heartbeat was detected. At the same time, the words "IH" (Irregular Heartbeat) will appear under the number of pulses recorded.

8. The measuring results will be shown until the **Start** button is pushed again or till automatic switch off 2 minutes later.

9. When the results was shown, press "M1" or "M2" key to record data in the Memory Set.

10. This device has average function. it will average the latest three measurements, which can help you to get your blood pressure values more correctly. Never rely on single blood pressure readings.

5.5 Fast interruption during a measurement


If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the **Start/Off** button can simply be pressed at any time. The device immediately reduces the cuff pressure.

5.6 Memory - displaying the last measurement


The device can store twin 60 sets (total of 120sets) of recorded data, and each set data includes values of diastolic and systolic pressure, pulse rate and date/time with lively figure. When "M1" is pressed, "01" will be displayed the symbols "M1" and "AVG" (average) and the average value of last three measurements recorded. To read the values of any measurement press again the M1 key, one time for any measurement (on the display will be displayed the record number of any measurement recorded, by more recently "1" to more old "2", "3", ... etc...)

Keep pressed for more than 2 seconds the M1 key to get data faster.


If your press the M1 key on he display will be displayed the "M1" symbol and its recorded data, if you press the "M2" key will be displayed the "M2" symbol and its recorded data.

If M1 is pressed If M2 is pressed

If M1 is pressed If M2 is pressed

NOTE: When the memories are full, the records will be erased one to one starting by more old when you record a new measurement. To erase all the recorded data you must take off the batteries for more than 10 seconds (or to disconnect the DC adapter by the unit).

Further information

Measurements should not occur soon after each other, since otherwise, the results will be falsified. Wait therefore for several minutes (at least 15 min.) in a relaxed position, sitting or lying, before repeating a measurement.

6. CARE, STORAGE AND REPAIR

- 1.It's necessary to protect this device against high moisture and heat, direct sunlight, shock, solvents, alcohol.
- 2.Remove the batteries if the device will be unused for a long time, and keep the batteries out of the reach of children.
- 3.Keep the cuff away from sharp subjects and don't extend or twist the cuff.
- 4.Use only soft and dry cloth to clean the device.
- 5.The cuff is delicate and must be handled with care. You can clean the cuff cover with damp cloth.

WARNING: Under no circumstances any part of device and its accessories could be washed.

Keep the device and all of its accessories out of the reach of children.

7. ERROR MESSAGE

| INDICATION | PROBLEM | CAUSE AND SOLUTION |
|---|---|--|
|  | Device doesn't measure. | Low Voltage or interior circuit problem. Remove and change the batteries. If the device is used with DC adapter, please check that the main connection plug is in compliance with rating indicated on the adapter. If after these checks the monitor still does not work, contact the customer service of Distributor. |
|  | Inflation not enough or deflation fault | <ol style="list-style-type: none"> 1. Cuff is not fastened well, please refasten it again according to the instruction. 2. Cuff is broken, change it with new one. 3. Pump, venting or deflation valve fail, please refer to the distributor. |
|  | Interferences and noises | <ol style="list-style-type: none"> 1. Move or talk during measurement. 2. Posture, etc. was different for each measurement. 3. User is too nervous. 4. Repeat the measurement following the manual and consult physician if the problem persists. |
| | Abnormal measurement | <p>When measuring, the value will increase and be shown on the display.</p> <ol style="list-style-type: none"> 1. If the value is not over 300 mmHg, check the good position of cuff and repeat the measure. If it still does not work please consult a physician. 2. If the value is over 300 mmHg, please contact the customer service of Distributor. |
| | Over safety pressure | <ol style="list-style-type: none"> 1. Air pipe blocked. 2. Inflation controller failure. 3. Pressure sensor fault. <p>Please check and repeat the measurement and/or ask to your distributor.</p> |
|  | Low Battery | Replace with four new batteries at the same time. |

8. TROUBLESHOOTING

| TROUBLE | CHECK POINT | SOLUTION |
|--|--|---|
| Though the batteries are installed, there is no indication on the display or the wrong indication. | <ol style="list-style-type: none"> 1. The batteries are not good positioned. 2. The batteries are exhausted. 3. The polarity of batteries is wrong. 4. The contacts of batteries compartment is dirty. | <ol style="list-style-type: none"> 1. Remove batteries and wait some minutes before install yet these. 2. Replace the batteries. 3. Install the batteries correctly. 4. Clean the battery terminals with dry cloth. |
| With the DC adapter installed, there is no indication on the display when the START button is pressed. | <ol style="list-style-type: none"> 1. The adapter is connected to a main socket with different rating. 2. The adapter is not well connected to the main socket or to the device. 3. The adapter is damaged. | <ol style="list-style-type: none"> 1. Check if data on the adapter are suitable for the main socket. 2. Check the connections of adapter to the device and main socket. 3. Insert batteries in the device and check if with batteries the device works. If yes, the adapter must be changed because damaged. |

| TROUBLE | CHECK POINT | SOLUTION |
|---|---|---|
| The inflation action cannot be done or the air pressure cannot rise. | 1. The cuff position is not correct. 2. The cuff connection is not correct. | 1. Check the cuff position and repeat the measurement. 2. Check the cuff connection and repeat the measurement. |
| The weak battery sign appears on the display. | Battery low. | Replace the batteries Important: In order to store the previous measurement data, ensure the replace action is done under the Clock mode. |
| The blood pressure cannot be taken and on the display appears error message. | 1. The cuff is not well fastened. 2. User is nervous. 3. The posture is not correct. | 1. Re-fasten the cuff. 2. Sit down and relax. Keep quiet during the measurement. 3. Keep the cuff and heart at the same height during the measurement. If the patient has severe heart problem, then the blood pressure may not be read correctly. |
| Under normal measurement circumstances, the reading at home is different from that of the clinics, or each measurement has various reading. | 1. The variation is due to the different environments. 2. The blood pressure changes with physiological or psychological status. | Record everyday value trend and consult it with a physician. |

9. TECHNICAL CHARACTERISTICS

| | |
|------------------------------------|--|
| Model | REALCHECK PRO 2 |
| Size | 131(W) x 128(H) x 86(D)mm |
| Weight | Approximately 400g not including batteries |
| Measuring method | Oscillometry |
| Measuring range | 20 to 280 mmHg (pressure) 40 to 200 beats/minute (pulse rate) |
| Measuring accuracy | ± 3 mmHg for systolic and diastolic pressure ± 5% of the reading for the pulse rate |
| Inflation | Automatic |
| Rapid deflation | Automatic valve |
| Data Transmission | Connection to PC with USB cable |
| Operating out | Automatic (after 2 minutes) |
| Batteries | 4x1.5V LR06 alkaline batteries type AA |
| Memory | Twin zones (60x2 memories) |
| Operation temperature and humidity | +10°C to + 40°C, 15% to 85%RH |
| Storage temperature and humidity | -5°C to + 50°C, 15% to 90%RH |
| Cuff size | Applicable for arm size 23-33 cm |

Notes: Specification are subject to change without further notice.

| Guidance and manufacture's declaration – Electromagnetic immunity (as request by regulation EN 60601-1-2:2001) | | | |
|--|------------------------------|------------------------------|---|
| The model Automatic Blood Pressure Monitor is tended for use in the electromagnetic environment specified below. The customer or the user of the model Automatic Blood Pressure Monitor should assure that is used in such an environment. | | | |
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environmet-guidance |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ± 6 kV contact ± 8 kV air | ± 6 kV contact ± 8 kV air | Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material the relative humidity should be at least 30%. |
| Electrical fast transient / burst IEC 61000-4-4 | Not applicable | - | - |
| Surge IEC 61000-4-5 | Not applicable | - | - |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | Not applicable | - | - |
| Power frequency magnetic field IEC 61000-4-8 | 3 A/m | 3 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment |

| Guidance and manufacture's declaration – Electromagnetic immunity (as request by regulation EN 60601-1-2:2001) | | | |
|--|---------------------------|------------------|---|
| The model REALCHEK PRO 2 Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of this device should assure that is used in such an environment. | | | |
| Immunity test | IEC / EN 60601 test level | Compliance level | Electromagnetic environmet-guidance |
| Conducted RF IEC 61000-4-6 | Not applicable | - | Portable and mobile RF communications equipment should be used no closer to any part of the REALCHECK PRO 2 Blood Pressure Monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. $d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80\text{MHz to } 800\text{MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800\text{MHz to } 2,5\text{GHz}$ <p>Where P is the maximum out put power rating of the transmitter in watts (W) according to the transmitter manufacturer and is recommended separation distance in meters (m). Field strength from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b). Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
| Ratiated RF IEC 61000-4-3 | 3V/m 80MHz-2.5GHz | 3V/m | |

Note 1 At 80MHz and 800MHz the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

(a) Field strength from fixed transmitters, such as base stations for radio telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the REALCHECK PRO 2 Blood Pressure Monitor should be observed to verify normal operation. If abnormal performance is observed additional measures may be necessary, such as re-orienting or relocating the device.

(b) Over the frequency range 150KHz to 80MHz, field strength should be less than 3V/m

10. SIMBOLOGY

| | |
|---|--|
|  | Type BF |
| CE 0123 | CE marking in conformity with EC directive 93/42/EEC |
|  | Warning, consult the instruction manual |
|  | To keep in dry and cool place |
|  | Storage temperature |
|  | Insulation Class II |
| ~ | Alternate Current |
| — — — | Direct Current |
| Start | Start |
| Off | Off |
| M1 - M2 | M1 Memory / Recall - M2 Memory / Recall |
| Mode | Mode of setting |
| Adjust | Data Adjust |
| Set | Setting |

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"Actuation of European directives 2002/95/EC, 2002/96/EC and 2003/108/EC, for reduction in use of dangerous substances in the electric and electronic device and for garbage disposal"

The symbol as over applied on the device or its packaging means that at the end of its useful life the product must not be disposed of with domestic waste. At the end of device useful, the user will must deliver it to the able collecting centres for electric and electronic garbage, or give back to the retailer in the moment of equivalent new device purchasing, one against one.

Disposing of the product separately prevents possible negative consequences for the environment and for health, deriving from inadequate disposal. It also allows the recovery of materials of which it's made up in order to obtain an important saving of energy and resources and to avoid negative effects to the ambient and health. In case of abusive disposal of device by user, will be applied administrative endorsements in compliance with current standard.

11. WARRANTY

Warranty conditions

- The unit is covered by a 2-year warranty against defects in materials and workmanship. The warranty does not extend to units that have undergone unauthorised repairs or been tampered with.
- The warranty period commences on the purchase date. It is mandatory to attach the sales receipt or invoice serving as proof of the purchase date to the warranty certificate below, which must be filled out in every part.
- Under the warranty, any components with manufacturing defects will be replaced and/or repaired free of charge. The Purchaser will bear all expenses for transporting, delivering and collecting the unit.

The warranty does not cover:

- The accessories supplied with the unit or parts subject to normal wear.
- Repair of unsubstantiated defects.
- Repairs of units that have failed or been rendered defective as a result of improper use, carelessness or negligence or if the damage is not attributable to the manufacturer (accidental falls, careless transport, etc.).

The warranty does not cover any direct or indirect damage or injury caused to persons, animals or property as a result of improper use of the product or during its period of inefficiency.

WHEN SUBMITTING A WARRANTY CLAIM, IN OBSERVANCE OF THE ABOVE CONDITIONS, SEND THE UNIT TOGETHER WITH THE CERTIFICATE BELOW AND PROOF OF PURCHASE TO:
SSL HEALTHCARE ITALIA S.P.A. VIA MARCO EMILIO LEPIDO, 178/5 - BOLOGNA (ITALY).

CERTIFICATE TO BE RETURNED IN CASE OF REPLACEMENT/REPAIRS

THE WARRANTY IS VALID ONLY IF THE SALES RECEIPT OR INVOICE IS ATTACHED

Unit Type/Model: _____

Lot: _____ Serial No: (SN) _____

Date of purchase: _____

PURCHASER'S DETAILS

Last Name and Name: _____

Address: _____

Telephone: _____

Description of the defect: _____

☐ I hereby authorise the use of the above information
in accordance with the Law on Privacy 196/03

Signature confirming acceptance of warranty conditions _____

— |

| —

— |

| —

CE 0123

Medical Device
Classe IIa
(DIR 93/42/CEE)

Fabbricante/Manufacturer:

CA.MI snc

Via Ugo La Malfa, 31
43010 Pilastro (PR) - Italy

Distributore/Distributor

SSL Healthcare Italia SPA

Via M.E.Lepido, 178/5
40132 Bologna - Italy